IMPORTANT
This information is intended to aid in using the Pioneer Sternal Cable System and is not a reference for surgical technique. Refer to the surgical technique manual for instructions for implantation. Only the instruments designed for the Pioneer Sternal Cable System should be used during the procedure.

SINGLE USE
All implants are single use and implanted system components should never be reused under any circumstances. An explanted device should never be re-implanted. If reused, single use devices may not perform as intended and could cause serious injury.

INDICATIONS
Cardiovascular sternal closure

CONTRAINDICATIONS
Contraindications may include, but are not limited to:

- Presence of documented infection;
- Patient metal allergy or intolerance;
- Presence of severe osteopenia and/or osteoporosis, rapid bone absorption, metabolic bone disease, cancer, tumor, or tumor like condition of the bone;
- Inadequate tissue coverage of implant site;
- Interference with other critical anatomical structures;
- Undiagnosed infection, end stage malignant disease, or other unexplained disease;
- Severely committed fractures;
- Any patient unwilling to follow postoperative instructions; and
- Any situation not defined by Indications.

WARNINGS AND PRECAUTIONS

- WARNING** Do not intermix implants with different metallic alloy types.
- WARNING** Do not use if packaging is damaged / opened prior to use.
- CAUTION** Storage damage due to poor understanding of the cable system, surgical technique, and bio-mechanical principles of cable fixation.
- CAUTION** Federal Law (USA) restricts the use of this device on the order of a physician only.
- IMPORTANT** Avoid over tensioning cables as they may break / fray.
- IMPORTANT** Cables may cut through soft bone that is not protected and immobilized.

This system has not been evaluated for safety and compatibility in the MR environment. This system has not been tested for heating or migration in the MR environment.

POSSIBLE ADVERSE EFFECTS
Irritation or inflammation of soft tissue structures surrounding implant

- Cables may cut through soft osteoporotic, osteopenic or cancellous bone if not properly protected and immobilized.
- Bone formation surrounding the implant (s) making removal difficult.
- Early fraying, kinking, loosening, or breaking of cable construct.
- Irritation from implants/ wires where there is not / inadequate soft tissue coverage over the implant.
- Loss of reduction of fractures or dislocations secondary to loosening, disassembly or breakage of the cables or secondary to cables cutting through osteoporotic bone.

Foreign body reaction causing possible tumor-like condition.

- Nonunion or delayed union of bone fracture or bony fusion.
- Rare possible neurovascular compromise, disruption of blood circulation, and/or storage damage due to improper cable placement and/or improper assembly of the systems components.
- Caution of growth of the operated portion of the bone.

NOTE: It may be necessary to perform additional surgery in order to correct some of the adverse effects or reactions that may or may not be related to the actual cable system.

PROOPERATIVE

- Avoid patients not meeting the criteria described above in the INDICATIONS portion of this package insert.
- Avoid patients with conditions that may predispose to a possible poor result or adverse effect.
- Use caution with damaged implants. Take care to prevent damage or environmental corrosion while handling, preparing, and storing components and instruments. Cable should not be unraveled, kinked or damaged and none of the components should be scratched prior to implantation.

NTRAOPERATIVE

- Instructions, manuals, protocols and/or videos should be carefully followed on the cable system surgical technique. Only the instruments designed for the cable system and the instruments used in the procedure should be used from the patient postoperatively.
- Cut off excess cable flush with the crimp at completion of the procedure and after fully crimping the crimp.
- Prior to the loosening of the tension in the tensioner-crimper device and prior to cutting the cable, care should be taken to make sure that the crimps are fully crimped according to the operative technique.

POSTOPERATIVE

It is extremely important that the physician provides clear directions and warnings and that the patient understands the conditions of the patient postoperatively:

- Partial- or non-weight bearing may be recommended or required to achieve firm bone union.
- Warn patient against sudden changes in position, strenuous activity, falls, smoking, consuming alcohol or other drugs not prescribed by the physician, steroids, non-steroidal anti-inflammatory agents, aspirin, and mechanical vibrations that may loosen the devices.
- If appropriate, restrict patient’s mobility at the fusion region.
- Follow all physician instructions carefully for a more successful and complete healing process.
- If the device is not removed after its functional purpose is complete, the following complications may occur:
 - Friction with loosened tissue or pain;
 - The implant may migrate leading to injury;
 - Injury may occur as result of postoperative trauma;
 - Breakage, loosening, fraying or make removal difficult;
 - Device presence may cause pain, discomfort, abnormal sensations, increased risk of infection; and/or bone loss due to stress shielding.

PRODUCT COMPLAINTS

Complaints or dissatisfaction with the cable system’s quality, safety, reliability, durability, effectiveness and/or performance, brought forth by a Health Care Professional, whether via a customer or user of the System of products, should be immediately conveyed to the attention of your distributor via telephone, FAX or written correspondence sent by express mail. It is important to note that when filing a complaint, the following information must be included in order to properly respond to the complaint:

Name and address; nature of the complaint; the component(s) trade name and catalog number; applicable lot number(s); and notification of whether or not a written report from the distributor is being solicited.

FURTHER INFORMATION
Recommended surgical operative techniques and uses for cable products are available at no charge upon request. Should subsequent information be requested or required, please contact the following:

INDICATIONS
Kardiiovaskulær sternal lukning

Pioneer Surgical Technology, Inc.
375 River Park Circle
Marquette, MI 49855 U.S.A.
Telefon: (906) 226-4457
Fax: (906) 226-4455

PHYSICIAN NOTE: General information that may be provided to the patient on the use and limitation of these devices is found below.

PATIENT INFORMATION

- Your implant is a temporary device designed to stabilize / secure your bone fracture(s) and augment the process of healing after which time, if conditions are unfavorable, your physician may choose to remove it.
- Avoid sudden changes in position, strenuous activity, falls, smoking, consuming alcohol or other drugs not prescribed by your physician, steroids, non-steroidal anti-inflammatory agents, aspirin, and mechanical vibrations that may loosen the devices.
- If you experience unusual pain, severe discomfort, or fever, contact your physician immediately.
- Follow all physician instructions carefully for a more successful and complete healing process.
- Cable products are comprised of titanium or stainless steel.

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PROOPERATIVE
Avoid patients not meeting the criteria described above in the INDICATIONS portion of this package insert.

- Avoid patients with conditions that may predispose to a possible poor result or adverse effect.
- Use caution with damaged implants. Take care to prevent damage or environmental corrosion while handling, preparing, and storing components and instruments. Cable should not be unraveled, kinked or damaged and none of the components should be scratched prior to implantation.

NTRAOPERATIVE
Instructions, manuals, protocols and/or videos should be carefully followed on the cable system surgical technique. Only the instruments designed for the cable system and the instruments used in the procedure should be used from the patient postoperatively.

Cut off excess cable flush with the crimp at completion of the procedure and after fully crimping the crimp.

Prior to the loosening of the tension in the tensioner-crimper device and prior to cutting the cable, care should be taken to make sure that the crimps are fully crimped according to the operative technique.

POSTOPERATIVE
It is extremely important that the physician provides clear directions and warnings and that the patient understands the conditions of the patient postoperatively:

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- Cables may cut through soft osteoporotic, osteopenic or cancellous bone if not properly protected and immobilized.
- Bone formation surrounding the implant (s) making removal difficult.
- Early fraying, kinking, loosening, or breaking of cable construct.
- Irritation from implants/ wires where there is not / adequate soft tissue coverage over the implant.
- Loss of reduction of fractures or dislocations secondary to loosening, disassembly or breakage of the cables or secondary to cables cutting through osteoporotic bone.

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