
## SINGLE WIRE & DOUBLEWIRE™ STERNAL SUTURES

### **PRODUCT DESCRIPTION:**

Sternal sutures are nonabsorbable, sterile, surgical sutures composed of 316L stainless steel. Stainless steel sutures are a permanent implant used in sternal closure, abdominal wound closure, hernia repair, and certain orthopedic procedures including cerclage and tendon repair. The suture consists of a length of monofilament wire and a mechanically attached surgical needle. The wire is implant grade 316L stainless steel and the needle is surgical grade stainless steel.

There are two configurations of Sternal Sutures, Single Wire and DoubleWire™. The sutures are available in multiple diameters, lengths, and configurations which vary according to surgical need and surgeon preference. A variety of needles of different sizes and tip configurations are attached to the wires also according to surgical need and surgeon preference.

### **INDICATIONS FOR USE:**

Stainless steel sutures are indicated for use by trained surgeons or surgical personnel in sternal closure, abdominal wound closure, hernia repair, and certain orthopedic procedures including cerclage and tendon repair. The basic operating principle is that of mechanical binding.

### **CONTRAINDICATIONS:**

The use of the suture is contraindicated in patients with known sensitivities or allergies to stainless steel, the primary constituents of which are chromium and nickel.

### **WARNINGS:**

1. As the risk of wound dehiscence may vary with the site of application and the suture material used, users should be familiar with surgical procedure and techniques involving nonabsorbable suture before using the product.
2. **DO NOT RE-USE SUTURES:** Known risks of re-use include device contamination, material degradation, malfunction, and patient/ user injury.

### **PRECAUTIONS:**

1. In handling the product, care should be taken to avoid damage from handling.
2. The sterile package should be inspected before use. Discard the package if the seal is open or the package is damaged.

### **POTENTIAL SIDE EFFECTS:**

Side effects associated with the use of this device include: wound dehiscence, suture breakage, allergic response in patients with known sensitivities to metals contained in 316L stainless steel (e.g. chromium and nickel), infected wounds, acute inflammatory tissue reaction, and transitory irritation at the wound site.

### **INSTRUCTIONS FOR USE:**

See precautions. Using aseptic technique, open the pouch, remove the holding card, open the card and remove the sutures. Refer to the following instructions for Single Wire and DoubleWire™ sutures:

#### Single Wire

1. Using the surgeon's preferred technique, the curved needle is used to pass the suture either para or trans sternally.
2. The curved needle is removed from the suture and the suture ends are twisted together using a needleholder and tightened to form a band around the sternum. Overtightening may stress the suture creating the potential for future breakage.
3. Using the surgeon's preferred technique, remove any excess twisted suture wire. Fold the remaining twisted end into the sternal seam or other location determined by the surgeon.


#### DoubleWire™

1. DoubleWire™ sutures are installed using the Wire Twister. For detailed instructions refer to the Wire Twister Instructions for Use.
2. When installing the DoubleWire™, the curved needle is passed around each side of the mesosternum so that an equal portion of the loop appears on both sides. Care should be taken to avoid the internal mammary artery. Single Wire sutures are placed through the manubrium and around the xiphoid. The curved needle is removed from the wire by cutting the single wire strand with wire cutters immediately above the DoubleWire™ sternal suture junction.
3. After installing the three DoubleWire™ sternum sutures and the two Single Wire sutures, the sternum is initially approximated with the single wires in the manubrium and xiphoid areas.
4. The hook of the DoubleWire™ Twister tool is passed through the two open portions of each loop on opposite sides of the sternum. Both loops must be located in the lowest portion of the twister hook prior to twisting.
5. Upward pressure is applied with the twister to draw each loop up evenly on each side of the sternum. While maintaining this upward pressure of the twister, then rotate to braid the two pair of wires together. The anterior plates of the sternum should be approximated with the index finger of the hand while twisting the wires.
6. The twisting tool is rotated until the two loops have been fully absorbed into a single cable.
7. At the onset of the double cabling, optimum tension has been achieved and twisting should be discontinued.
8. The wires are cut with double action wire cutters (i.e. Aesculap #LX156) leaving about ¼ to 1 cm.
9. Then turn down the twisted wire in the direction of the twist to bury the tip of the wire. The twisted wire is then flattened onto the anterior plates of the sternum with the blunt tip of a needleholder.

**STORAGE:**

The sterile packages should be stored at room temperature not to exceed 50°C (120° F).

PRODUCT LABEL SYMBOLS			
	Do not use if package is damaged		Reverse Cutting Needle Tip
	Single Suture		Conventional Cutting Needle Tip
	Double Suture		Taper Cutting Needle Tip
	1/2 Circle Needle		Diamond Point Needle Tip
	3/8 Circle Needle		Ball End Needle Tip

The Instructions for Use and a comprehensive symbol glossary can be found at <http://www.aemedical.com/quality/>

	<b>A&amp;E MEDICAL CORPORATION</b> 5206 Asbury Road, P.O. Box 758 Farmingdale, NJ 07727 USA	<b>EC</b> <b>REP</b>	<b>CE</b> 0344
		BIOMED, S.A. Parque Tecnológico de Madrid C/ Einstein, 3 Tres Cantos 28760 Madrid, Spain	

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