PRODUCT DESCRIPTION:

MYOWire® Ultra-Thin Temporary Pacing Wire is an electrical connection between the heart and equipment used to monitor the heart's activity or provide external pacing signals to the heart. The electrical connection consists of a positive and a negative lead. The leads are attached to the heart surface during surgical procedures and the opposite ends are passed through the chest wall exterior to the body. The exterior end is then connected to the monitoring/pacing equipment using an appropriate patient cable. Based on patient need and user preference, one of the leads may be placed externally on the skin to complete the electrical connection. Once the electrical connection is complete, monitoring and pacing is possible.

MYOWire® Ultra-Thin Temporary Pacing Wires are offered in various configurations to address patient and physician needs and preferences as well as facilitating different attachment methods. Atrial attachments (epicardial) are configured for surface suturing. Ventricular attachments (intramyocardial) are configured for intramuscular insertion with various curved needle options. Skin attachments (indifferent leads) are configured with either a straight or curved needle for inserting through the skin.

INDICATIONS FOR USE:

MYOWire® Ultra-Thin Temporary Pacing Wires are intended for temporary atrial and ventricular pacing and sensing for contemplated implant duration of 7 days or less. The device is supplied sterile, non-pyrogenic and intended for single use.

CONTRAINDICATIONS:

There are no known general contraindications to temporary peri- or post-surgical pacing. The particular medical condition and anatomy of the patient, however, may dictate the wire system and implantation procedure to be used. Temporary pacing wires should not be used when an internal pacemaker is implanted in the patient.

PRECAUTIONS:

1. Ensure that hospital personnel involved in the care and use of temporary pacing wires are aware of the warnings and precautions associated with temporary pacing wires.
2. Standard hospital procedures concerning maintaining sterility and preventing infection are to be followed.
3. Defibrillation equipment should be readily available during installation and removal of the temporary pacing wire.
4. Removal of the pacing wire should only be performed by trained, experienced healthcare professionals who are aware of the potential risks associated with removal.

WARNINGS:

1. An inserted pacing wire forms a direct current path to the myocardium. Therefore, avoid the possibility of unintentional contact between the pacing wires and any equipment. During handling (insertion, connecting, testing) of pacing wires, the terminal pins or exposed metal are not to be touched nor be allowed to contact electrically conductive or wet surfaces, to protect against fibrillation that may be caused by alternating currents. During pacing wire insertion and testing, use only battery powered equipment (or line powered equipment specially designed for this purpose).
2. Line powered equipment used in the vicinity of the patient must be properly grounded. Pacing wire connectors must be insulated from any leakage currents that may arise from line powered equipment.
3. Do not reuse. Known risks of re-use include device contamination, material degradation, malfunction, and patient/user injury.
4. The MYOWire® Ultra-Thin Temporary Pacing Wires are for single use only. Do not re-use or re-sterilize the product.
5. Use of a pacing wire in a Magnetic Resonance environment (Magnetic Resonance Imaging examination) is not recommended.
6. While transporting a patient, the pacing wire needs to be secured to avoid the possibility of the pacing wire becoming dislodged.

POTENTIAL COMPLICATIONS:

The possibility of dislodgement, wire fractures, threshold elevation and sensing inconsistency exists. Other potential complications include, but are not limited to, myocardial irritability, septicemia and myocardial infections, myocardial (heart tissue) bleeding and damage, especially during removal of the pacing wire.

INSTRUCTIONS FOR USE:

Remove the thoracic needle and enough wire from the spool to pass the needle through the abdominal wall. Pull the wire through the abdominal wall and remove all the wire from the spool leaving enough wire within the thorax cavity to attach the electrodes to the heart. Remove the heart needle from the spool and attach the electrodes to the heart.


The physician may choose to vary the placement and fixation procedure in accordance with his best judgment. The selected fixation method should always allow for later pacing wire removal which is to be performed by applying gentle traction. The product shall only be used by trained physicians or staff appointed by a physician.

FIXATION METHODS:

Run the heart needle through the myocardium (Fig. 1) in a direction that the entry point of the pacing wire in the heart can later be aligned with the exit from the chest. This eases the withdrawal of the wire when removing it.

- The distal electrode and anchor flares are pulled through the myocardial tissue and exposed at the site of the needle pathway (Fig. 2). The distal electrode is then pulled back into the myocardial tissue which will open out the anchor flares onto the epicardial surface.
- It is important that both electrodes (distal and proximal) are fully in the myocardium.

- Cut off the heart needle leaving a minimum of 1 cm but not more than 2 cm of wire to the distal electrode (Fig. 2, 3, 4).
- The surgeon may choose not to remove the heart needle immediately as there is a risk of dislodgement.
- The heart needle has to be removed before surgery is completed.
- Should the wire become dislodged after the curved heart needle has been cut off but before closure of the chest, repositioning should not be attempted. A new wire should be used.
- Before closing the patient, the pacing wire should be tested with an external pacemaker to verify signal transfer to the heart tissue.
- If a pacing wire fails twice to operate after verification of proper placement and connection to the external pacemaker, replace the wire with a new sterile wire.

If no heart needle is attached, the stripped end of the pacing wire should be sutured directly onto the myocardium (Fig. 5, 6).

**CONNECTION TO EQUIPMENT**

Hold the thoracic needle (distal to the breakaway marking) and run it through the chest wall in such a direction that the pacing wire in the heart is aligned with the exit from the chest. Leave an extra length of the pacing wire in the chest to prevent dislodgement during patient movement.


1. Snap the thoracic needle at the breakaway marking and discard the needle.
2. Slide the thoracic needle from the contact pins and discard the needle.
3. Verify pacing wire to external pulse generator compatibility before use
4. Connect the pacing wire to the appropriate terminals on the external pulse generator.

REMOVAL

When the pacing wire is no longer required prior to hospital discharge (typically within 7 days of surgery), gentle traction should be applied to remove the wire. Precautions should be considered in accordance with the judgment and experience of the physician.

SELF STRIPPING ADAPTOR

1. Each strand of the pacing wire requires an individual adaptor.
2. Slide each strand of the wire into the hole of plastic part A (Fig. 7) leaving a minimum of 2 cm wire in the adaptor. All pieces of the strand shall be kept in the adaptor.
3. Close the adaptor by folding the adaptor halves together until the two locks snap in place.
4. The protective cap on the metal adaptor pins can be used to prevent unintentional contact.

**ADAPTOR USED WITH CONTACT PINS**

Each strand of the pacing wire requires an individual adaptor. Place the end of the metal heart wire connector pin inside the lumen of the adaptor pin (Fig. 8). Ensure the heart wire connector pin is completely inside the plastic adaptor body, at exact position of the bridge (see red arrow Fig. 8). If not positioned on the bridge there is no conductive contact, before closing the adaptor lid. Close the adaptor lid by pressing the adaptor lid. Ensure it locks into the main body of the adaptor.

To re-open adaptor with the lid of the adaptor housing facing upwards (Fig.9), hold the adaptor connector pin with one hand while using the thumb of the other hand to press the raised part at the rear of the adaptor housing backwards and down. This will allow the adaptor lid to un-lock. It can be opened, and the heart wire connector pin removed. Warning: if heart wire connector pin is not positioned on the bridge there is no conductive contact.

**STORAGE:**

The sterile package should be stored at room temperature.

PRODUCT LABEL SYMBOLS	
	Number of needles in pouch.

The Instructions for Use and a comprehensive symbol glossary can be found at <http://www.aemedical.com/quality/>

	A&E MEDICAL CORPORATION 5206 Asbury Road, P.O. Box 758 Farmingdale, NJ 07727 USA	EC	REP	
		BIOMED, S.A. Parque Tecnológico de Madrid C/ Einstein, 3 Tres Cantos 28760 Madrid, Spain		

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