EN

## ROTATING SURGICAL PUNCH

**PRODUCT DESCRIPTION:**

The Rotating Surgical Punch is a sterile, single use, disposable, surgical punch which creates a circular hole in soft tissue. The Rotating Surgical Punch consists of a 3 piece (body, finger grip, & thumb rest) plastic housing shaped to fit the user's hand. The circular cut is created by the instrument grade stainless steel die and pin set. Various sized cutting diameters between 2.7 mm and 5.0 mm are available for the Rotating Surgical Punch to enable it to be used according to surgical need.

**INDICATIONS FOR USE:**

The surgical punch is intended for use by trained surgeons to create circular cuts in soft tissue or openings in soft vessel walls during surgical procedures. The basic operating principle is that of mechanical shearing using round dies and blades.

**CONTRAINDICATIONS:**


None known

**WARNINGS:**

1. Only surgeons trained in the particular procedure should use the punch.
2. After activation, do not release the plunger until the punch is removed from the incision site to prevent the excised tissue from being ejected into the surgical area or vessel interior.
3. **DO NOT RE-USE.** Known risks of re-use include device contamination, material degradation, malfunction, and patient/user injury.

**INSTRUCTIONS FOR USE:**

1. Inspect the sterile package. Do not use if the package is opened or damaged. Open the package using sterile, operating room procedure.
2. Remove the protective cap from the punch tip. Examine the punch for damage paying particular attention to the tip and blade. Do not use if damage is present.
3. The punch is activated by placing your fingers under the punch handles (B) and your thumb on the plunger (A). Pressing your thumb and fingers together causes the punch tip (E) to move towards the punch body (C) past the circular cutting blade (D).
4. Activate the punch several times prior to surgical use. Do not use if the punch does not move freely.
5. Make an incision in the area where the punch is to be used. The incision should be slightly smaller than the punch tip. Insert and center the punch tip through the incision.
6. Holding the punch perpendicular to the tissue being cut, activate the punch. Do not release the plunger.
7. Keeping the plunger depressed, remove the punch from the incision area.
8. Slowly release the plunger to ensure that the excised tissue remains captured beneath the punch tip.
9. Remove the excised tissue from the punch tip.

**STORAGE:**

The sterile package should be stored at room temperature not to exceed 50°C (120°F).

PRODUCT LABEL SYMBOLS			
	Open pouch in this direction		Do not use if package is damaged

The Instructions for Use and a comprehensive symbol glossary can be found at <http://www.aemedical.com/quality/>

	<b>A&amp;E MEDICAL CORPORATION</b> 5206 Asbury Road, P.O. Box 758 Farmingdale, NJ 07727 USA	<table border="1"> <tr> <td>EC</td> <td>REP</td> </tr> </table>	EC	REP	0344
		EC	REP		
BIOMED, S.A. Parque Tecnológico de Madrid C/ Einstein, 3 Tres Cantos 28760 Madrid, Spain					

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