


MYO/KLIP™ Connecting System

Catalog
Number




Description: **MYO/KLIP™** connector protects exposed distal end of temporary pacing wire.

MYO/KLIP™ Sterile, 12/Box
MYO/KLIP™ Non-Sterile, 120/Box

019-510
019-511


MYO/LEAD™ Cable



MYO/LEAD™ cable accepts **MYO/KLIP™** connector at one end and plugs into temporary pacing unit at other end. (Medtronic 5320, 5342, 5345, 5375, 5880, 5880A)

MYO/LEAD™ 6' (1.8m)

019-501



MYO/LEAD™ cable accepts **MYO/KLIP™** connector at one end and plugs into temporary pacing unit at other end. (Medtronic 5348 and 5388)

019-503

*Non-stock product. Available as custom order product in minimum 10 box quantities.