







# MYO/WIRE™ Temporary Cardiac Pacing Wires

Indifferent Leads	Length	Insulation Color	Reverse Cutting 	Taper Point 
 <p data-bbox="62 559 942 678">Description: Eight inch ground or indifferent electrode. Special curved cutting needle permits subcutaneous placement of electrode.</p>			 <p data-bbox="1605 616 1750 688">28mm 1/2 Circle</p>	 <p data-bbox="1854 616 1989 688">41mm Straight</p>
 <p data-bbox="62 932 735 968">Electrode 1.0mm Diameter x 22mm Length</p>	<p data-bbox="1108 714 1274 740">8" (20cm)</p> <p data-bbox="1108 756 1274 782">8" (20cm)</p> <p data-bbox="1108 797 1274 823">8" (20cm)</p> <p data-bbox="1108 859 1274 885">24" (61cm)</p> <p data-bbox="1108 901 1274 927">24" (61cm)</p> <p data-bbox="1108 942 1274 968">24" (61cm)</p>	<p data-bbox="1336 714 1460 740">Orange</p> <p data-bbox="1336 756 1440 782">White</p> <p data-bbox="1336 797 1512 823">Lt. Orange</p> <p data-bbox="1336 859 1460 885">Orange</p> <p data-bbox="1336 901 1440 927">White</p> <p data-bbox="1336 942 1512 968">Lt. Orange</p>	<p data-bbox="1605 714 1740 740">027-523</p> <p data-bbox="1605 756 1740 782">027-524</p> <p data-bbox="1605 797 1740 823">027-525</p> <p data-bbox="1605 859 1740 885">027-520</p> <p data-bbox="1605 901 1740 927">027-521</p> <p data-bbox="1605 942 1740 968">027-522</p>	<p data-bbox="1854 714 1989 740">027-043</p> <p data-bbox="1854 756 1989 782">027-044</p> <p data-bbox="1854 797 1989 823">027-045</p> <p data-bbox="1895 859 1947 885">----</p> <p data-bbox="1895 901 1947 927">----</p> <p data-bbox="1895 942 1947 968">----</p>